

INVERTED decompression

DFM™

DECOMPRESSION & FUNCTIONAL MOVEMENT
Inversion Table



Experience the versatility of the DFM™ - Decompression & Functional Movement Inversion Table.

Maximizing function in a small space, the DFM™ is an effective therapeutic tool for any clinic. Patients may rest prone or supine while inverted, allowing for treatments at varied angles. Adjustable features ensure user comfort while the easily accessible push-button control allows for supervised or unassisted inversion sessions.

Teeter™

DFM™

DECOMPRESSION & FUNCTIONAL MOVEMENT
Inversion Table

The **slanted foam bed** provides an effective stretching platform and is designed with a **detachable foam piece** to allow for comfortable treatment in the prone position.

Adjustable knee bolsters provide knee flexion when supine and support the legs when prone.

Adjustable forearm rests can be adapted to fit to the height of any user and help to achieve proper positioning during treatment.

The knee-high **easy-reach ankle closure** opens with the simple push of an ergonomically shaped handle, assisting patients with limited mobility. To ensure user security, the **weighted handle design** closes automatically upon table rotation should this vital step be overlooked.

The **gravity lock ankle clamp system**, constructed of aerospace-grade stainless steel parts, serves as a secondary mechanism to prevent accidental opening of the ankle clamps in the inverted position.

A **durable UL-Listed motor** with push-button control enables smooth and accurate rotation to any angle. **Interlocking gears** in the motor enable bed stabilization at any angle for hassle-free treatment.



DIMENSIONS: 35" x 46" x 61" (89 x 117 x 155 cm)
WEIGHT: 180 lbs. (81.8 kg)
INCLUDED: Training Manual
WARRANTY: Two Years

Teeter